



Ideal for Gingival Scalloping & Sinus Window Preparation

SURGICAL SCIENCE KERRISON RONGEUR

*Special Price**
\$295

reg.: \$395

CALL 1.877.668.0888

* Offer expires December 31, 2007



Kerrison Rongeur Scalloping Gingiva



1. Edentulous upper lateral at 2nd stage.



2. Reverse-action cutting end of KR



3. Crestal incision with KR under buccal gingival flap.



4. Scalloped gingival contour



5. Prepared gingiva contacting temporary abutment

Sinus Window preparation



1. Cutting buccal sinus wall with KR



2. Serial bone cuts to enlarge window

Features of the Surgical Science Kerrison Rongeur:

- Delicate reverse action bone cutting for sinus window opening.
- Low profile head for safe access and visibility in the sinus.
- Creation of attached gingival scalloped margin around implant abutments.
- 7.6cm shaft.
- 4mm cut.
- This high quality instrument has a **10 year** full replacement warranty on breakage!